

# Bee9: An M9 derived medium for bee gut bacteria

## Description

This protocol is used to prepare stock solutions for an M9-based 'minimal' medium containing trace metals and vitamins, without casamino acids. Additional vitamins, and growth nutrients are added for growth of *S. alvi* and/or *Gilliamella* strains. The pH should be adjusted to 5.5

---

### 1. Preparation of HMB (hutner's mineral base) modified metals 44 (component 1)

- Weight on a balance the following:
  1. 1.095g  $\text{ZnSO}_4 \times 7\text{H}_2\text{O}$ , dissolve in **5ml** ddH<sub>2</sub>O at pH=2
  2. 0.914g  $\text{FeSO}_4 \times 7\text{H}_2\text{O}$ , dissolve in **5ml** ddH<sub>2</sub>O at pH=2
  3. 0.154g  $\text{MnSO}_4 \times \text{H}_2\text{O}$ , dissolve in **5ml** ddH<sub>2</sub>O at pH=2
  4. 0.392g  $\text{CuSO}_4 \times 5\text{H}_2\text{O}$ , dissolve in **100ml** ddH<sub>2</sub>O at pH=2
  5. 0.248g  $\text{Co}(\text{NO}_3)_2 \times 6\text{H}_2\text{O}$ , dissolve in **100ml** ddH<sub>2</sub>O at pH=2
  6. 0.177g  $\text{Na}_2\text{B}_4\text{O}_7 \times 10\text{H}_2\text{O}$ , dissolve in **100ml** ddH<sub>2</sub>O at pH=2
- Mix **1.**, **2.**, **3.**, and add 10ml of each of **4.**, **5.**, **6.** to obtain **component 1.**
- Adjust volume to 100ml and sterile-filter with 0.22 um filter
- Keep it at 4C, light protected

### 2. Preparation of additional salts (components 2, 3, 4, 5)

- **Component 2:** 1.72g of  $\text{MgSO}_4 \times 7\text{H}_2\text{O}$ , dissolve in **100ml** ddH<sub>2</sub>O at pH=2, sterile-filter with 0.22 um filter
- **Component 3:** 3.33g of  $\text{CaCl}_2 \times 2\text{H}_2\text{O}$ , dissolve in **100ml** ddH<sub>2</sub>O at pH=2, sterile-filter with 0.22 um filter
- **Component 4:** 0.99g of  $\text{FeSO}_4 \times 7\text{H}_2\text{O}$ , dissolve in **100ml** ddH<sub>2</sub>O at pH=2, sterile-filter with 0.22 um filter
- **Component 5:** 0.974g of  $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \times 4\text{H}_2\text{O}$ , dissolve in **100ml** ddH<sub>2</sub>O at pH=2, sterile-filter with 0.22 um filter
- Keep at RT, light protected

### 3. Preparation of the vitamin 10x stock solution (component 6a)

- 0.2g of Calcium pantothenate, dissolve in **10ml** ddH<sub>2</sub>O
- 0.01g of Thiamine-HCl, dissolve in **10ml** ddH<sub>2</sub>O
- 0.01g of biotin, dissolve in **10ml** ddH<sub>2</sub>O
- Mix 79ml of ddH<sub>2</sub>O, 10ml of **2.**, 10ml of **3.**, and 1ml of **4.**
- Sterile-filter the final 100 ml of vitamin stock solution (component 6a)
- Component 6a can be kept for at least 6 months at 4°C in the fridge, protected from light

### 5. Preparation of 10x stock solution of M9 (component 8)

- Mix the following:
  - 60g  $\text{Na}_2\text{HPO}_4$
  - 30g  $\text{KH}_2\text{PO}_4$
  - 5g  $\text{NaCl}$
  - 10g  $\text{NH}_4\text{Cl}$
- Add water to 1l, autoclave
- Keep at RT

## 6. Preparation of *S. alvi* specific vitamins

- Weigh out powder of the following compounds and dissolve in M9 base to a concentration of 10 mM (for final concentration of [0.1 mM])
- 4-Hydroxybenzoate
- P-aminobenzoic acid
- Pyridoxine\*HCl
- Choline

## 7. Prepare 500ml of M9 base medium

- Prepare the final M9 base medium freshly on the day of use
- For 500ml mix the following:
- Component 1: 500 µl
- Component 2: 8390 µl
- Component 3: 1000 µl
- Component 4: 100 µl
- Component 5: 10 µl (100 uL of 1/10 component 5 for small M9 volume)
- Component 6a (10x Vitamin stock): 50ml
- Component 8 (10x M9 stock): 50ml
- Custom vitamins (10mM stock) 5 mL
- ddH<sub>2</sub>O: 385 ml

## 8. Prepare carbon + media

- Dissolve compounds of choice in M9 to a final concentration of 10 mM. Citrate is a good single carbon source. Glucose [10 mM] can also be added for *Gilliamella* growth.
- Common substrates include.
  - Alpha-ketoglutaric acid
  - Fumaric acid
  - Malic acid
  - Succinic acid
  - Citric acid
  - Sodium pyruvate
  - Sodium acetate
  - Sodium lactate
- Check pH and adjust to 5.5 to 6

